




WHITE PAPER

Magnetic Module for Bead-Based Extraction and Purification

Magnetic Plate to Engage and Disengage with Labware

Written by
Opentrons



SECTION 1

Product Description

The Opentrons Magnetic Module is a device that automatically engages and disengages high-strength neodymium magnets to seated well plates within the module housing. The module can readily be used with the Opentrons App. The Magnetic Module comes with two separate adjustable plate brackets for supporting standard and deep well labware. Users are able to remove and/or swap brackets by unscrewing the module's thumb screw.

The magnetic plate on the module supports many different magnetic bead products for extraction and purification of nucleic acids in the labware shown below. The brackets allow for support of protocols with standard or deep well 96 plates with volumes ranging from 10 μ L to 2 mL. For more information in designing a protocol on the Magnetic Module, please refer to the [Nucleic Acid Purification Technical Note](#).



Magnetic Module with 0.2 mL 96 well PCR plate




Magnetic Module with 2 mL Deep Well Plate

SECTION 2

Magnetic Module Overview

WELL PLATE COMPATIBILITY

The standard plate bracket supports 0.2 mL plate sizes up to 22 mm in vertical height. The deep well bracket can support deep well trays up to 47 mm in vertical height.



SECTION 3

Data

EVALUATION IN NEST FULL SKIRT 96 WELL PLATES

The Opentrons GEN2 Magnetic Module has been evaluated to determine the amount of time it takes for the supernatant to become clear. All trials took place at the default height (18mm) on the magnetic module for 8 bead types. Each test was performed in a NEST full skirt 96 well plate. For each bead type, 4 volumes were measured; 20 μ l, 30 μ l, 40 μ l and 50 μ l. At each volume, wells A1, A6, A12, D6, E7, H1 and H12 were filled to determine how pelleting changed across the plate. Starting from 0 minutes (Figure 1), images of the plate were taken in 1 minute increments and the wells were observed until all wells had a clear supernatant (Figure 2). The results for the extraction clean up beads can be seen in Table 1, the results for the NGS kit beads can be seen in Table 2 and the results for library preparation can be seen in Table 3.

FIGURE 1

96 well plate with no pelleting




FIGURE 2

96 well plate with clear supernatant



EVALUATION IN NEST DEEP WELL PLATES

For the deep well plate experiments, 1 ml and 2 ml deep well plates were used. The same designated wells were used as before, in the 96 well plates on the Opentrons GEN2 Magnetic Module at the default height for each test. For these experiments, only Omega-Biotek's Mag-Bind HDQ beads were used. To prepare these beads, 400 μ L HDQ Binding Buffer and 20 μ L Mag-Bind® Particles HDQ were added to each well and pipetted repeatedly to mix well. As with the 96 well plate, starting from 0 minutes (Figure 3), images of the plate were taken in 1 minute increments and the wells were observed until all wells had a clear supernatant.

FIGURE 3

Nest 1 mL Deep well plate with beads fully suspended




FIGURE 4

Nest 1 mL Deep well plate with clear supernatant




TABLE 1

Pelleting Times for Extraction Clean Up Beads at Default Height

| COMPANY | BEADS | VOLUME | ENGAGE HEIGHT | TIME TO PELLET |
|-----------------|-----------------------|--------|---------------|----------------|
| Beckman Coulter | SPRISelect | 20 µL | default | 6 min |
| Beckman Coulter | SPRISelect | 30 µL | default | 6 min |
| Beckman Coulter | SPRISelect | 40 µL | default | 8 min |
| Beckman Coulter | SPRISelect | 50 µL | default | 8 min |
| Qiagen | Qiaseq | 20 µL | default | 5 min |
| Qiagen | Qiaseq | 30 µL | default | 5 min |
| Qiagen | Qiaseq | 40 µL | default | 5 min |
| Qiagen | Qiaseq | 50 µL | default | 7 min |
| Omega Bio-tek | Mag-bind RXNPure Plus | 20 µL | default | 5 min |
| Omega Bio-tek | Mag-bind RXNPure Plus | 30 µL | default | 5 min |
| Omega Bio-tek | Mag-bind RXNPure Plus | 40 µL | default | 6 min |
| Omega Bio-tek | Mag-bind RXNPure Plus | 50 µL | default | 8 min |

TABLE 2

Pelleting Times for NGS Kit Beads at Default Height

| COMPANY | BEADS | VOLUME | ENGAGE HEIGHT | TIME TO PELLET |
|-----------------|------------------------|--------|---------------|----------------|
| Beckman Coulter | AMPure XP | 20 µL | default | 5 min |
| Beckman Coulter | AMPure XP | 30 µL | default | 5 min |
| Beckman Coulter | AMPure XP | 40 µL | default | 5 min |
| Beckman Coulter | AMPure XP | 50 µL | default | 6 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 20 µL | default | 5 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 30 µL | default | 5 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 40 µL | default | 6 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 50 µL | default | 8 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 20 µL | default | 5 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 30 µL | default | 5 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 40 µL | default | 6 min |
| Omega Bio-tek | Mag-Bind TotalPure NGS | 50 µL | default | 8 min |
| BASE Prime | PureBASE NGS | 20 µL | default | 6 min |
| BASE Prime | PureBASE NGS | 30 µL | default | 6 min |
| BASE Prime | PureBASE NGS | 40 µL | default | 7 min |
| BASE Prime | PureBASE NGS | 50 µL | default | 7 min |

TABLE 3

Pelleting Times for Library Preparation Beads at Default Height

| COMPANY | BEADS | VOLUME | ENGAGE HEIGHT | TIME TO PELLET |
|----------|--|--------|---------------|----------------|
| iGenomix | Riptide High Throughput Rapid Library Prep | 20 µL | default | 4 min |
| iGenomix | Riptide High Throughput Rapid Library Prep | 30 µL | default | 6 min |
| iGenomix | Riptide High Throughput Rapid Library Prep | 40 µL | default | 7 min |
| iGenomix | Riptide High Throughput Rapid Library Prep | 50 µL | default | 8 min |
| illumina | Nextera DNA Flex Library Prep | 20 µL | default | 4 min |
| illumina | Nextera DNA Flex Library Prep | 30 µL | default | 4 min |
| illumina | Nextera DNA Flex Library Prep | 40 µL | default | 7 min |
| illumina | Nextera DNA Flex Library Prep | 50 µL | default | 6 min |

TABLE 4

Deep Well Results

| MAGNET | PLATE | COMPANY | BEADS | VOLUME | ENGAGE HEIGHT | TIME TO PELLET |
|--------|----------------|---------------|--------------|--------|---------------|----------------|
| N42 | 1 mL Deep Well | Omega Bio-tek | Mag-Bind HDQ | 420 µL | default | 1 min |
| N42 | 2 mL Deep Well | Omega Bio-tek | Mag-Bind HDQ | 420 µL | default | 1 min |

CONSIDERATIONS

When using heavy beads such as Omega-Biotek's Mag-Bind HDQ beads, it is important to note that these beads are heavy and settle quickly when left to stand. Beads should be mixed thoroughly immediately prior to usage to prevent premature settling of the beads from the supernatant.

POWER USAGE

Module input: 36V, 2A

Power adapter input: 100-240VAC, 50/60Hz

SECTION 4


Dimensional Drawings

Standalone Magnetic Module Mass is 1.5kg.

Millimeter units are shown in all dimensional drawings.


DIMENSIONAL DRAWINGS

Magnetic Module



DIMENSIONAL DRAWINGS

Magnetic Module with Standard Plate Bracket



DIMENSIONAL DRAWINGS

Magnetic Module with Deep Well Plate Bracket

