Mikro-Dispenser, fix
Micro dispenser, fix
Micro-distributeur, fixe
Microdosificador, fijo

Mikro-Dispenser, umschaltbar
Micro dispenser, dual volume
Micro-distributeur, á 2 vol.
Microdosificador, dual

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Gebrauchsanleitung
Operating Manual
Mode d'emploi
Instrucciones de manejo

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Before using the instrument for the first time, ensure it is rinsed carefully or discard the first few samples dispensed.

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Contents

1. Safety Instructions	4
2. Application and Operating Limitations	5
3. Operating Exclusions	6
4. Operating Limitations	7
5. Components	8
6. Dispensing	9
6.1. Preparation for dispensing	9
6.2. Priming	9
6.3. Adjusting the volume	10
6.4. Dispensing	10
7. Cleaning	11
7.1. Emptying	11
7.2. Standard cleaning	12
7.3. Intensive cleaning	13
8. Cleaning/replacing the valves	14
8.1. Intake valve	14
8.2. Discharge valve	15
9. Calibration	16
10. Troubleshooting	18
11. Repair service	19
12. Warranty	19
13. Technical Data/Ordering Information	20
13.1. Items supplied	20
13.2. Accessories and Spare Parts	21

Table des matières		Contenido	
1. Règles de sécurité	22	1. Normas de seguridad	22
2. Fonction et restrictions d'emploi	23	2. Función y limitaciones de uso	23
3. Interdictions d'ensemble	24	3. Excepciones de uso	24
4. Restrictions d'emploi	25	4. Limitaciones de uso	25
5. Dessin d'ensemble	26	5. Dibujo de conjunto	26
6. Dosage	27	6. Dosificación	27
6.1. Préparatifs pour le dosage	27	6.1. Preparar la dosificación	27
6.2. Désaéragé	27	6.2. Purgar el aire	27
6.3. Changer le volume	28	6.3. Cambiar el volumen	28
6.4. Dosage	28	6.4. Dosificar	28
7. Nettoyage	29	7. Limpieza	29
7.1. Vidange	29	7.1. Vaciar	29
7.2. Nettoyage standard	30	7.2. Limpieza estándar	30
7.3. Nettoyage intensif	31	7.3. Limpieza a fondo	31
8. Nettoyage/remplacement des soupapes	32	8. Limpieza/cambio de las vlvulas	32
8.1. La soupape d'aspiration	32	8.1. La vlvula de aspiracin	32
8.2. La soupape d'jection	33	8.2. La vlvula de expulsin	33
9. Calibrage	34	9. Calibracin	34
10. Drangement	35	10. Averas	36
11. Retour pour rparation	37	11. Envos para reparacin	37
12. Garantie	37	12. Garanta	37
13. Donnees techniques/ Donnees de commande	38	13. Datos tcnicos/ Referencias	38
13.1. Emballage standard	38	13.1. Alcance del suministro	38
13.2. Accessoires et pices de rechange	39	13.2. Accesorios y piezas de recambio	39

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1. Safety Instructions

This Manual does not purport to address every safety issue which may arise during use. It is the responsibility of whomever uses this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Use the utmost caution when dispensing caustic, poisonous, radioactive or hazardous chemicals.
2. Observe general safety regulations (e.g., wear protective clothing, goggles and gloves).
3. Observe the Operating Manual and information from reagent manufacturers.
4. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations. If in doubt regarding the suitability of the instrument, for a particular application contact the manufacturer.
5. Before use, always verify that the instrument is in good working order, e.g., piston moves smoothly; filling and discharge tubes are firmly seated and are properly attached, etc.
6. Never use force on the instrument. Use of force may result in injury to the user or other persons.

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2. Application and Operating Limitations

This instrument is designed for dispensing liquids, observing the following physical limits:

- 15 to 40 °C for instrument and reagent
- density up to 2.2 g/cm³
- vapor pressure up to 500 mbar
- viscosity up to 200 mm²/s

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3. Operating Exclusions

When the instrument is correctly used, the dispensed liquid comes into contact with only the following materials:

Borosilicate glass, FEP, ETFE, PFA, PTFE, Platinum-Iridium and PVDF (closure cap).

The instrument is suitable for dispensing a wide range of liquids, **excluding**:

- Liquids which attack FEP, ETFE, PFA, PTFE and PVDF;
- Solutions containing hydrofluoric acid, since they attack borosilicate glass;
- Suspensions, containing solid particles;
- Solutions which decompose and form solid particles (e.g., Biuret reagent);
- Substances which undergo catalytic transformation or react with Platinum-Iridium (e.g., H_2O_2);
- Carbon disulfide

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Liquids which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or highly concentrated alkaline solutions).

When dispensing flammable media, make provisions to avoid static charging, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

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Warning!


If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions (see chapter 7) or contact the manufacturer.

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Note:

Compatibility of the instrument for this special application (e.g., trace material analysis) must be checked by the user or contact the manufacturer.

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6. Dispensing

6.1. Preparation for dispensing

1. Observe Operating Limitations and general safety regulations.
2. Push the intake tube into the intake valve as far as possible and cut the lower end at an angle.
3. Screw the dispenser on to the reagent bottle and align it with the bottle label.

Warning!

Avoid splashing of reagent!

4. Place a suitable vessel below the discharge tube.
5. The discharge tube must point away from the user at all times.

6.2. Priming

1. Place a suitable vessel below the discharge tube.

Warning!

Reagent may splash!


2. Gently push down the piston several times until air bubbles have disappeared from the discharge tube.

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6.3. Adjusting the volume (only piccolo 2)

Warning!

Reagent will dispense while volume is being adjusted. Hold a suitable vessel below the discharge tube outlet.

To set the volume, push down the operating knob and twist handle 90°. The nose of the housing indicates the volume after releasing operating knob.

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6.4. Dispensing

1. Place a suitable vessel below the discharge tube.
2. Slowly steadily depress piston, release gently. No air bubbles should be visible during dispensing.

Warning!

If the piston moves stiffly or is difficult to move, stop dispensing immediately and clean the instrument (see chap. 7).

Attention!

Reagent may escape from the discharge tube.

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7. Cleaning

Instrument will only function safely if cleaned regularly. Be sure to clean the instrument:

1. Immediately, if the piston moves stiffly;
2. Before changing reagents
3. Before longer breaks in use
4. Before any maintenance and repair work

Warning!

Dispenser, intake tube and discharge tube may be filled with reagent. To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Always point the discharge tube outlet away from the user or other people.

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7.1. Emptying

1. While the dispenser is still mounted on the bottle, place into a suitable basin.
2. Unscrew the dispenser and lift it high enough so that the intake tube is no longer immersed in liquid.
3. Cautiously tap the intake tube against the inside of the bottle so that any remaining reagent runs out.
4. Remove the dispenser from the bottle and mount it on another empty bottle.

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5. Hold the discharge tube outlet over the opening of the first bottle. Empty the discharge tube by repeated dispensing movements.

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7.2. Standard cleaning

1. Mount the empty dispenser on a bottle filled with a suitable cleaning solution.

Warning!

Reagent may splash!

2. Rinse the dispenser by repeated dispensing movements. Place a suitable vessel below the discharge tube.
3. Completely empty the instrument by repeated dispensing movements, as described in chapter 7.1.
4. Mount the dispenser on a bottle filled with distilled water and rinse it by repeated dispensing movements. Then empty the instrument as described in chapter 7.1.

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7.3. Intensive cleaning

Intensive cleaning must be preceded by a standard cleaning. This procedure must be followed if the piston is difficult to move. For this purpose, the dispenser must be partially disassembled.

Warning!

Before disassembling, always perform the standard cleaning procedure. To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Avoid splashing of reagent.

1. Pull off the intake tube (11) and clean it with a soft bottle brush. Replace if necessary, especially if cracked or worn.
2. Use a small screw driver to gently pry between the securing ring slot (3) and valve head (8) and carefully push up the securing ring (3).
3. Disassemble outer sleeve (2) by pressing down operating knob (4) and pulling the securing ring (3) simultaneously.
4. Remove the operating knob (4) with the piston (5) and the "Lifting spring".
5. Clean the piston and the cylinder with a soft bottle brush and rinse with water.

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Attention!

Insert the piston at a slight angle to avoid damage to the lip seal.

6. Mount the lifting spring (6) and carefully reinsert the piston (5).
7. Replace the outer sleeve (2) and press it down until it audibly locks in place.
8. Push down the securing ring (3).

8. Cleaning/replacing the valves

8.1. Cleaning/replacing the intake valve

Note:

Valve ball (9) may fall out when valve is removed

1. Unscrew the intake valve (10) using the valve tool (17).
2. Clean or replace the valve.
3. Make sure that valve ball is still in place. To reassemble, handtighten valve (10) and complete using valve tool (17).

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8.2. Cleaning/replacing the discharge valve

1. Remove the discharge tube support (12).
Pull out securing clip (13) with a small screwdriver.
2. Pull out the discharge valve (14).
3. Unscrew the discharge tube securing nut (15) and pull out the discharge tube (16).
4. Clean or replace the valve.
5. To reassemble push in valve (14) by hand and insert securing clip (13).
6. Push in the discharge tube as far as possible and tighten the securing nut (15).
7. In reverse order mount the discharge tube support (12).

Warning!

Never use force during disassembly and assembly. Verify that all components fit tightly and securely. Subsequent to reassembly, gravimetrically check volume (see chapter 9).

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9. Calibration

The procedure is described in detail in DIN EN ISO 8655/6

1. Set instrument to nominal volume and dispense distilled water into a weighing vessel.
2. Weigh the dispensed amount of water with an analytical balance.
3. Convert the indicated weight into volume units. taking into account the temperature, density and atmospheric pressure. (The conversion factor Z is published in the standard.)
4. Carry out 10 dispensing procedures.
5. Calculate systematic error (A%) and random error (CV%).

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10. Troubleshooting

Trouble	Possible causes	Action to be taken
Piston is difficult to move	Formation of crystals	Clean instrument (see chapter 7)
Reagent is not aspirated	Sticking filling valve	Clean intake valve; replace if necessary (see chapter 8.1)
Cylinder does not fill	Discharge valve damaged	Clean discharge valve; replace if necessary (see chapter 8.2)
Intake-/discharge tube not firmly connected	Tubes damaged	Replace intake-/discharge tube
Air is aspirated	Intake valve not firmly connected	Tighten intake valve; replace if necessary (see chapter 8.1)
	Intake tube not firmly connected or damaged	Push intake tube on firmly; replace if necessary
	Reagent with high vapor pressure has been drawn in too quickly	Aspirate reagent more slowly
Dispensed volume is too low	Intake valve leaks	Clean intake valve; replace if necessary (see chapter 8.1)
	Intake tube not firmly connected or damaged	Push intake tube on firmly; replace if necessary

If these measures do not help, return instrument for repair (see chap. 11).

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11. Repair Service

Attention!

Only cleaned and decontaminated instruments can be evaluated/ repaired.

Clean and decontaminate the instrument thoroughly. Return instrument along with a description of the malfunction and dispensed liquids.

Complete the 'Declaration on Absence of Health Hazards' (ask your supplier or manufacturer for the form. The form can also be downloaded from www.vitlab.com).

Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used. Return shipment is carried out at the risk and expense of the sender.

12. Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

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13. Technical data / Ordering information

13.1 Items supplied

VITLAB® piccolo 1 or
VITLAB® piccolo 2 complete with
discharge tube, intake tube and tool.

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Number of testing procedures:

10 according to
DIN EN ISO 8655

(Technical specifications subject
to change!)

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**13.2 Accessories and Spare Parts
for VITLAB® piccolo 1 + 2**

Intake tube length 150 mm	1650010
Intake valve complete with ball (ETFE/Borosilicate glass)	1655011
Discharge tube	1650145
Discharge valve (PFA/Borosilicate glass/ Platinum-Iridium)	1655090
Threaded bottles, 100 ml, coated, GL28	1671505

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