


Brilliance in a Dispenser Syringe



The Choice for Dispenser and Pump Syringes

- Exceptional performance
- Longer life
- Precise fit
- Greater reliability
- Compatible with existing systems
- Custom design options

www.sge.com


SGE Dispenser and Pump Syringes incorporate Diamond Technology, taking advantage of new and innovative technologies and bringing those benefits to your laboratory.

For over 50 years, SGE syringes have been designed and manufactured for a variety of instruments and customer applications.

At the heart of your pump or dispenser, and combined with the rest of the system components, syringes must provide the superior performance you demand and rely upon.

Designed and optimized with outstanding quality and precise fit, SGE syringes provide greater reliability for applications run on your systems.

The Choice for Dispenser and Pump Syringes





SGE Diamond Pump and Dispenser Syringes:

- Improve the accuracy and precision of your results
- Are robust and durable for longer life, less maintenance downtime and reduced maintenance costs
- Do not contaminate the aspirated fluid
- Are designed to perfectly fit your instrument making them easy to use and install
- Increase flexibility with a broad functional temperature range and high chemical resistance
- Reduce downtime due to leakage
- Are available with 30 mm and 60 mm stroke lengths and capacities from 12.5 μ L to 100 mL
- Have an inert flow path. Only high quality borosilicate glass and premium PTFE contact aspirated fluids



Look for the diamond mark - the mark of syringe brilliance.



Type	Capacity	Description	Part Number
30 mm stroke, 1/4-28 UNF thread connection, 1/4" hole in plunger – For use with V3 pumps			
	50 µL	50D-V3-GT Dispenser Syringe	004693
	100 µL	100D-V3-GT Dispenser Syringe	005693
	250 µL	250D-V3-GT Dispenser Syringe	006693
	500 µL	500D-V3-GT Dispenser Syringe	007693
	1.0 mL	1MD-V3-GT Dispenser Syringe	008193
	2.5 mL	2.5MD-V3-GT Dispenser Syringe	008593
	5.0 mL	5MD-V3-GT Dispenser Syringe	008793
60 mm stroke, 1/4-28 UNF thread connection, 1/4" hole in plunger – For use with V6 pumps			
	25 µL	25D-V6-GT Dispenser Syringe	003696
	50 µL	50D-V6-GT Dispenser Syringe	004696
	100 µL	100D-V6-GT Dispenser Syringe	005696
	250 µL	250D-V6-GT Dispenser Syringe	006696
	500 µL	500D-V6-GT Dispenser Syringe	007696
	1.0 mL	1MD-V6-GT Dispenser Syringe	008196
	2.5 mL	2.5MD-V6-GT Dispenser Syringe	008596
	5.0 mL	5MD-V6-GT Dispenser Syringe	008796
	10 mL	10MD-V6-GT Dispenser Syringe	008996
	25 mL	25MD-V6-GT Dispenser Syringe	009496
30 mm stroke, Luer Lock connection – For use with PSD/3 and PSD/4 pumps			
	12.5 µL	12.5D-TLLX30-GT Dispenser Syringe	002692
	25 µL	25D-TLLX30-GT Dispenser Syringe	003694
	50 µL	50D-TLLX30-GT Dispenser Syringe	004694
	100 µL	100D-TLLX30-GT Dispenser Syringe	005694
	125 µL	125D-TLLX30-GT Dispenser Syringe	005692
	250 µL	250D-TLLX30-GT Dispenser Syringe	006694
	500 µL	500D-TLLX30-GT Dispenser Syringe	007694
	1 mL	1MD-TLLX30-GT Dispenser Syringe	008194
	1.25 mL	1.25MD-TLL30-GT Dispenser Syringe	008192
	2.5 mL	2.5MD-TLL30-GT Dispenser Syringe	008594
	5.0 mL	5MD-TLL30-GT Dispenser Syringe	008794
	12.5 mL	12.5MD-TLL30-GT Dispenser Syringe	008994
60 mm stroke, 1/4-28 UNF thread connection, 6-32 UNC thread in plunger button, CX has a button stop			
	50 µL	50D-CX-GT 50UL Syringe	004995
	100 µL	100D-CX-GT 100UL Syringe	005990
	250 µL	250D-CX-GT 250UL Syringe	006995
	500 µL	500D-CX-GT 500UL Syringe	007995
	1.0 mL	1MD-C-GT 1ML Syringe	008185
	2.5 mL	2.5MD-C-GT 2.5ML Syringe	008687
	5.0 mL	5MD-C-GT 5ML Syringe	008787
	10 mL	10MD-C-GT 10ML Syringe	008987
60 mm stroke, Luer Lock connection, 6-32 UNC thread in plunger button			
	1.0 mL	1MDF-LL-GT 1ML Syringe	008025
	2.5 mL	2.5MDF-LL-GT 2.5ML Syringe	008425
	5.0 mL	5MDF-LL-GT 5ML Syringe	008762
	5.0 mL	5MDR-LL-GT 5ML Syringe	008760
	10 mL	10MDF-LL-GT 10ML Syringe	008962
	10 mL	10MDR-LL-GT 10ML Syringe	008960
	25 mL	25MDF-LL-GT 25ML Syringe	009463
	25 mL	25MDR-LL-GT 25ML Syringe	009462
	50 mL	50MR-LL-GT 50ML Syringe	009660
	100 mL	100MR-LL-GT 100ML Syringe	009760

For more information visit www.sge.com or contact techsupport@sge.com.