




designed for scientists



## MX carrier

/// Data Sheet

The Matrix plate attachment is suitable for microtiter and deep-well plates. The ideal attachment design ensures an optimal temperature transfer and thus guarantees valid results. The thermoshakers Matrix Orbital Delta plus and Matrix Delta plus automatically recognize, which attachment is connected to the device via an interface in the locking mechanism. The attachment can be changed quickly and easily with an ergonomic lever.

Speed/ temperature max.:

3000 rpm (microtiter plates); 100 °C

2000 rpm (deep-well plates);

# IKA

designed for scientists